

PCT

WORLD INTELLECTUAL PROPERTY ORGANIZATION  
International Bureau



INTERNATIONAL APPLICATION PUBLISHED UNDER THE PATENT COOPERATION TREATY (PCT)

(51) International Patent Classification<sup>7</sup> :

A61F 2/06

A1

(11) International Publication Number:

WO 00/07523

(43) International Publication Date:

17 February 2000 (17.02.00)

(21) International Application Number: PCT/CA99/00695

(22) International Filing Date: 29 July 1999 (29.07.99)

(30) Priority Data:

60/094,950

31 July 1998 (31.07.98)

US

(71) Applicant (for all designated States except US): NOVO RPS  
ULC [CA/CA]; 865 West 10th Avenue, Vancouver, British  
Columbia V5Z 1L7 (CA).

(72) Inventors; and

(75) Inventors/Applicants (for US only): RICCI, Donald, R.  
[CA/CA]; 4443 West 3rd Avenue, Vancouver, British  
Columbia V6R 1M9 (CA). SHUKOV, George, A. [US/US];  
14440 De Bell Road, Los Altos Hills, CA 94022 (US).  
PENN, Ian, M. [CA/CA]; 6360 Larch Street, Vancouver,  
British Columbia V6R 4E9 (CA).

(74) Agents: NASSIF, Omar, A. et al.; Gowling, Strathy &  
Henderson, Suite 4900, Commerce Court West, Toronto,  
Ontario M5L 1J3 (CA).

(81) Designated States: AE, AL, AM, AT, AU, AZ, BA, BB, BG,  
BR, BY, CA, CH, CN, CU, CZ, DE, DK, EE, ES, FI, GB,  
GD, GE, GH, GM, HR, HU, ID, IL, IN, IS, JP, KE, KG,  
KP, KR, KZ, LC, LK, LR, LS, LT, LU, LV, MD, MG, MK,  
MN, MW, MX, NO, NZ, PL, PT, RO, RU, SD, SE, SG, SI,  
SK, SL, TJ, TM, TR, TT, UA, UG, US, UZ, VN, YU, ZA,  
ZW, ARIPO patent (GH, GM, KE, LS, MW, SD, SL, SZ,  
UG, ZW), Eurasian patent (AM, AZ, BY, KG, KZ, MD,  
RU, TJ, TM), European patent (AT, BE, CH, CY, DE, DK,  
ES, FI, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE), OAPI  
patent (BF, BJ, CF, CG, CI, CM, GA, GN, GW, ML, MR,  
NE, SN, TD, TG).

Published

With international search report.

(54) Title: BIFURCATED STENT DELIVERY SYSTEM AND METHOD OF USE

(57) Abstract

An endovascular sleeve which can be utilized to navigate a pair of guidewires to a bifurcated body passageway such that, once in place, the guidewires are substantially untwisted or untangled. This greatly facilitates delivery of the bifurcated stent to the bifurcated artery.

